

Galectin-9 Polyclonal Antibody

Catalog No: #40947



Package Size: #40947-1 50ul #40947-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Galectin-9 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	Galectin-9 Polyclonal Antibody detects endogenous levels of Galectin-9 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Galectin-9.
Target Name	Galectin-9
Other Names	LGALS9; Galectin-9; Gal-9; Ecalectin; Tumor antigen HOM-HD-21
Accession No.	Swiss-Prot: O00182NCBI Gene ID: 3965
Uniprot	O00182
GeneID	3965;
SDS-PAGE MW	40kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.


ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of HepG2 cells using Galectin-9 Polyclonal Antibody



Western Blot analysis of VEC cells using Galectin-9 Polyclonal Antibody

Note: This product is for in vitro research use only