

Gasdermin C Polyclonal Antibody

Catalog No: #40951

Package Size: #40951-1 50ul #40951-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Gasdermin C Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	Gasdermin C Polyclonal Antibody detects endogenous levels of Gasdermin C protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Gasdermin C.
Target Name	Gasdermin C
Other Names	GSDMC; MLZE; Gasdermin-C; Melanoma-derived leucine zipper-containing extranuclear factor
Accession No.	Swiss-Prot: Q9BYG8NCBI Gene ID: 56169
Uniprot	Q9BYG8
GeneID	56169;
SDS-PAGE MW	55kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of 293 cells using Gasdermin C Polyclonal Antibody

Note: This product is for in vitro research use only