

GATA-4 Polyclonal Antibody

Catalog No: #40954

Package Size: #40954-1 50ul #40954-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GATA-4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	GATA-4 Polyclonal Antibody detects endogenous levels of GATA-4 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human GATA-4 around the non-phosphorylation site of S105.
Target Name	GATA-4
Other Names	GATA4; Transcription factor GATA-4; GATA-binding factor 4
Accession No.	Swiss-Prot: P43694NCBI Gene ID: 2626
Uniprot	P43694
GeneID	2626;
SDS-PAGE MW	48kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HeLa cells using GATA-4 Polyclonal Antibody

Note: This product is for in vitro research use only