

GPR119 Polyclonal Antibody

Catalog No: #40973

Package Size: #40973-1 50ul #40973-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GPR119 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF ELISA
Species Reactivity	Hu
Specificity	GPR119 Polyclonal Antibody detects endogenous levels of GPR119 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human GPR119.
Target Name	GPR119
Other Names	GPR119; Glucose-dependent insulinotropic receptor; G-protein coupled receptor 119
Accession No.	Swiss-Prot: Q8TDV5NCBI Gene ID: 139760
Uniprot	Q8TDV5
GeneID	139760;
SDS-PAGE MW	37kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of COLO cells using GPR119 Polyclonal Antibody

Note: This product is for in vitro research use only