

GPR126 Polyclonal Antibody

Catalog No: #40975

Package Size: #40975-1 50ul #40975-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GPR126 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu
Specificity	GPR126 Polyclonal Antibody detects endogenous levels of GPR126 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human GPR126.
Target Name	GPR126
Other Names	GPR126; DREG; VIGR; G-protein coupled receptor 126; Developmentally regulated G-protein-coupled receptor; Vascular inducible G protein-coupled receptor
Accession No.	Swiss-Prot: Q86SQ4NCBI Gene ID: 57211
Uniprot	Q86SQ4
GeneID	57211;
SDS-PAGE MW	145kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/5000.

Not yet tested in other applications.

Images



Western Blot analysis of HuvEc cells using GPR126 Polyclonal Antibody

Note: This product is for in vitro research use only