

GPR152 Polyclonal Antibody

Catalog No: #40976

Package Size: #40976-1 50ul #40976-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	GPR152 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF
Species Reactivity	Hu
Specificity	GPR152 Polyclonal Antibody detects endogenous levels of GPR152 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human GPR152.
Target Name	GPR152
Other Names	GPR152; PGR5; Probable G-protein coupled receptor 152; G-protein coupled receptor PGR5
Accession No.	Swiss-Prot: Q8TDT2NCBI Gene ID: 390212
Uniprot	Q8TDT2
GeneID	390212;
SDS-PAGE MW	51kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunofluorescence: 1/200 - 1/1000.

Images



Western Blot analysis of k562 cells using GPR152 Polyclonal Antibody



Western Blot analysis of various cells using GPR152 Polyclonal Antibody

Note: This product is for in vitro research use only