

GPR158 Polyclonal Antibody

Catalog No: #40977



Package Size: #40977-1 50ul #40977-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GPR158 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms
Specificity	GPR158 Polyclonal Antibody detects endogenous levels of GPR158 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human GPR158.
Target Name	GPR158
Other Names	GPR158; KIAA1136; Probable G-protein coupled receptor 158
Accession No.	Swiss-Prot: Q5T848NCBI Gene ID: 57512
Uniprot	Q5T848
GeneID	57512;
SDS-PAGE MW	150kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/5000.

Not yet tested in other applications.

Images



Western Blot analysis of COLO 293T HELA cells using GPR158 Polyclonal Antibody

Note: This product is for in vitro research use only