

GPR173 Polyclonal Antibody

Catalog No: #40979

Package Size: #40979-1 50ul #40979-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GPR173 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	GPR173 Polyclonal Antibody detects endogenous levels of GPR173 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human GPR173.
Target Name	GPR173
Other Names	GPR173; SREB3; Probable G-protein coupled receptor 173; Super conserved receptor expressed in brain 3
Accession No.	Swiss-Prot: Q9NS66NCBI Gene ID: 54328
Uniprot	Q9NS66
GeneID	54328;
SDS-PAGE MW	47kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of 3T3 cells using GPR173 Polyclonal Antibody

Note: This product is for in vitro research use only