

GPR175 Polyclonal Antibody

Catalog No: #40980

Package Size: #40980-1 50ul #40980-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	GPR175 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF
Species Reactivity	Hu Ms Rt
Specificity	GPR175 Polyclonal Antibody detects endogenous levels of GPR175 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human GPR175.
Target Name	GPR175
Other Names	TPRA1; GPR175; PP6566; Transmembrane protein adipocyte-associated 1; Integral membrane protein GPR175
Accession No.	Swiss-Prot: Q86W33NCBI Gene ID: 131601
Uniprot	Q86W33
GeneID	131601;
SDS-PAGE MW	41kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Immunofluorescence: 1/200 - 1/1000.

Images



Western Blot analysis of 293 cells using GPR175 Polyclonal Antibody

Note: This product is for in vitro research use only