

GPR40 Polyclonal Antibody

Catalog No: #40981

Package Size: #40981-1 50ul #40981-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GPR40 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF ELISA
Species Reactivity	Hu Mk
Specificity	GPR40 Polyclonal Antibody detects endogenous levels of GPR40 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human GPR40.
Target Name	GPR40
Other Names	FFAR1; GPR40; Free fatty acid receptor 1; G-protein coupled receptor 40
Accession No.	Swiss-Prot: O14842NCBI Gene ID: 2864
Uniprot	O14842
GeneID	2864;
SDS-PAGE MW	26kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunofluorescence: 1/200 - 1/1000.


ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of COS7 cells using GPR40 Polyclonal Antibody



Western Blot analysis of COS-7 cells using GPR40 Polyclonal Antibody

Note: This product is for in vitro research use only