

Histone H3.3 Polyclonal Antibody

Catalog No: #41019

Package Size: #41019-1 50ul #41019-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Histone H3.3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	Histone H3.3 Polyclonal Antibody detects endogenous levels of Histone H3.3 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Histone H3.3 around the non-phosphorylation site of S31.
Target Name	Histone H3.3
Other Names	H3F3A; H3.3A; H3F3; PP781; H3F3B; H3.3B; Histone H3.3
Accession No.	Swiss-Prot: P84243NCBI Gene ID: 3020/3021
Uniprot	P84243
GeneID	3020;3021;
SDS-PAGE MW	15kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1:500-2000

Immunohistochemistry: 1/100 - 1/300

Immunofluorescence: 1/200 - 1/1000

ELISA: 1/5000

Images



Western Blot analysis of Hela cells using Histone H3.3 Polyclonal Antibody

Note: This product is for in vitro research use only