

HNF-4 α Polyclonal Antibody

Catalog No: #41023



Package Size: #41023-1 50ul #41023-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HNF-4 α Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	HNF-4 α Polyclonal Antibody detects endogenous levels of HNF-4 α protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human HNF4- α around the non-phosphorylation site of S313.
Target Name	HNF-4 α
Other Names	HNF4A; HNF4; NR2A1; TCF14; Hepatocyte nuclear factor 4-alpha; HNF-4-alpha; Nuclear receptor subfamily 2 group A member 1; Transcription factor 14; TCF-14; Transcription factor HNF-4
Accession No.	Swiss-Prot: P41235NCBI Gene ID: 3172
Uniprot	P41235
GeneID	3172;
SDS-PAGE MW	52kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of HT29 cells using HNF-4α Polyclonal Antibody

Note: This product is for in vitro research use only