

HoxA6 Polyclonal Antibody

Catalog No: #41029

Package Size: #41029-1 50ul #41029-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HoxA6 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	HoxA6 Polyclonal Antibody detects endogenous levels of HoxA6 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human HoxA6.
Target Name	HoxA6
Other Names	HOXA6; HOX1B; Homeobox protein Hox-A6; Homeobox protein Hox-1B
Accession No.	Swiss-Prot: P31267NCBI Gene ID: 3203
Uniprot	P31267
GeneID	3203;
SDS-PAGE MW	26kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of K562 cells using HoxA6 Polyclonal Antibody

Note: This product is for in vitro research use only