

Hrs Polyclonal Antibody

Catalog No: #41035

Package Size: #41035-1 50ul #41035-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hrs Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	Hrs Polyclonal Antibody detects endogenous levels of Hrs protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Hrs.
Target Name	Hrs
Other Names	HGS; HRS; Hepatocyte growth factor-regulated tyrosine kinase substrate; Hrs; Protein pp110
Accession No.	Swiss-Prot: O14964NCBI Gene ID: 9146
Uniprot	O14964
GeneID	9146;
SDS-PAGE MW	86kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of Jurkat cells using Hrs Polyclonal Antibody

Note: This product is for in vitro research use only