

HSP90β Polyclonal Antibody

Catalog No: #41041

Package Size: #41041-1 50ul #41041-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HSP90β Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IP IF ELISA
Species Reactivity	Hu Ms Rt Mk
Specificity	HSP90β Polyclonal Antibody detects endogenous levels of HSP90β protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human HSP90β around the non-phosphorylation site of S254.
Target Name	HSP90β
Other Names	HSP90AB1; HSP90B; HSPC2; HSPCB; Heat shock protein HSP 90-beta; HSP 90; Heat shock 84 kDa; HSP 84; HSP84
Accession No.	Swiss-Prot: P08238NCBI Gene ID: 3326
Uniprot	P08238
GeneID	3326;
SDS-PAGE MW	83kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300. Immunoprec


Immunoprecipitation: 2-5 ug/mg lysate.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of NIH-3T3 cells using HSP90 β Polyclonal Antibody

Note: This product is for in vitro research use only