

HXK I Polyclonal Antibody

Catalog No: #41049

Package Size: #41049-1 50ul #41049-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HXK I Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	HXK I Polyclonal Antibody detects endogenous levels of HXK I protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human HXK I.
Target Name	HXK I
Other Names	HK1; Hexokinase-1; Brain form hexokinase; Hexokinase type I; HK I
Accession No.	Swiss-Prot: P19367NCBI Gene ID: 3098
Uniprot	P19367
GeneID	3098;
SDS-PAGE MW	109kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.


ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HeLa cells using HXK I Polyclonal Antibody



Western Blot analysis of MOUSE-BRAIN cells using HXK I Polyclonal Antibody

Note: This product is for in vitro research use only