

Integrin β5 Polyclonal Antibody

Catalog No: #41071



Package Size: #41071-1 50ul #41071-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Integrin β5 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt Mk
Specificity	Integrin β5 Polyclonal Antibody detects endogenous levels of Integrin β5 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Integrin β5.
Target Name	Integrin β5
Other Names	ITGB5; Integrin beta-5
Accession No.	Swiss-Prot: P18084NCBI Gene ID: 3693
Uniprot	P18084
GeneID	3693;
SDS-PAGE MW	88kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

Images



Western Blot analysis of COS7 cells using Integrin β5 Polyclonal Antibody



Western Blot analysis of Jurkat cells using Integrin $\beta 5$
Polyclonal Antibody

Note: This product is for in vitro research use only