

JDP2 Polyclonal Antibody

Catalog No: #41083

Package Size: #41083-1 50ul #41083-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	JDP2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	JDP2 Polyclonal Antibody detects endogenous levels of JDP2 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human JDP2 around the non-phosphorylation site of T148.
Target Name	JDP2
Other Names	JDP2; Jun dimerization protein 2
Accession No.	Swiss-Prot: Q8WYK2NCBI Gene ID: 122953
Uniprot	Q8WYK2
GeneID	122953;
SDS-PAGE MW	22kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of Jurkat cells using JDP2 Polyclonal Antibody

Note: This product is for in vitro research use only