

KDEL Receptor 2 Polyclonal Antibody

Catalog No: #41090

Package Size: #41090-1 50ul #41090-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	KDEL Receptor 2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	KDEL Receptor 2 Polyclonal Antibody detects endogenous levels of KDEL Receptor 2 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human KDEL Receptor 2.
Target Name	KDEL Receptor 2
Other Names	KDEL2; ERD2.2; ER lumen protein retaining receptor 2; ERD2-like protein 1; ELP-1; KDEL endoplasmic reticulum protein retention receptor 2; KDEL receptor 2
Accession No.	Swiss-Prot: P33947NCBI Gene ID: 11014
Uniprot	P33947
GeneID	11014;
SDS-PAGE MW	24kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of COLO205 cells using KDEL Receptor 2 Polyclonal Antibody

Note: This product is for in vitro research use only