

KIR2.1 Polyclonal Antibody

Catalog No: #41093

Package Size: #41093-1 50ul #41093-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	KIR2.1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Rt
Specificity	KIR2.1 Polyclonal Antibody detects endogenous levels of KIR2.1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human KIR2.1.
Target Name	KIR2.1
Other Names	KCNJ2; IRK1; Inward rectifier potassium channel 2; Cardiac inward rectifier potassium channel; Inward rectifier K(+) channel Kir2.1; IRK-1; hIRK1; Potassium channel; inwardly rectifying subfamily J member 2
Accession No.	Swiss-Prot: P63252NCBI Gene ID: 3759
Uniprot	P63252
GeneID	3759;
SDS-PAGE MW	48kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of HeLa cells using KIR2.1 Polyclonal Antibody

Note: This product is for in vitro research use only