

LGR6 Polyclonal Antibody

Catalog No: #41106

Package Size: #41106-1 50ul #41106-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	LGR6 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF ELISA
Species Reactivity	Hu
Specificity	LGR6 Polyclonal Antibody detects endogenous levels of LGR6 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human LGR6.
Target Name	LGR6
Other Names	LGR6; Leucine-rich repeat-containing G-protein coupled receptor 6
Accession No.	Swiss-Prot: Q9HBX8NCBI Gene ID: 59352
Uniprot	Q9HBX8
GeneID	59352;
SDS-PAGE MW	104kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of HuvEc cells using LGR6 Polyclonal Antibody

Note: This product is for in vitro research use only