

MDMX Polyclonal Antibody

Catalog No: #41128

Package Size: #41128-1 50ul #41128-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MDMX Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	MDMX Polyclonal Antibody detects endogenous levels of MDMX protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human MDMX around the non-phosphorylation site of S367.
Target Name	MDMX
Other Names	MDM4; MDMX; Protein Mdm4; Double minute 4 protein; Mdm2-like p53-binding protein; Protein Mdmx; p53-binding protein Mdm4
Accession No.	Swiss-Prot: O15151NCBI Gene ID: 4194
Uniprot	O15151
GeneID	4194;
SDS-PAGE MW	45kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of 293T cells using MDMX Polyclonal Antibody

Note: This product is for in vitro research use only