

MEF-2 Polyclonal Antibody

Catalog No: #41130

Package Size: #41130-1 50ul #41130-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MEF-2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	MEF-2 Polyclonal Antibody detects endogenous levels of MEF-2 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human MEF-2 around the non-phosphorylation site of S408.
Target Name	MEF-2
Other Names	MEF2A; MEF2; Myocyte-specific enhancer factor 2A; Serum response factor-like protein 1
Accession No.	Swiss-Prot: Q02078NCBI Gene ID: 4205
Uniprot	Q02078
GeneID	4205;
SDS-PAGE MW	55kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of HeLa cells using MEF-2 Polyclonal Antibody

Note: This product is for in vitro research use only