

MLKL Polyclonal Antibody

Catalog No: #41150



Package Size: #41150-1 50ul #41150-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	MLKL Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	MLKL Polyclonal Antibody detects endogenous levels of MLKL protein.
Immunogen Description	Synthesized peptide derived from the N-terminal region of human MLKL.
Target Name	MLKL
Other Names	MLKL; Mixed lineage kinase domain-like protein
Accession No.	Swiss-Prot: Q8NB16NCBI Gene ID: 197259
Uniprot	Q8NB16
GeneID	197259;
SDS-PAGE MW	55kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Images



Western Blot analysis of K562 cells using MLKL Polyclonal Antibody

Note: This product is for in vitro research use only