

Mnk1 Polyclonal Antibody

Catalog No: #41160

Package Size: #41160-1 50ul #41160-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mnk1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms
Specificity	Mnk1 Polyclonal Antibody detects endogenous levels of Mnk1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Mnk1.
Target Name	Mnk1
Other Names	MKNK1; MNK1; MAP kinase-interacting serine/threonine-protein kinase 1; MAP kinase signal-integrating kinase 1; MAPK signal-integrating kinase 1; Mnk1
Accession No.	Swiss-Prot: Q9BUB5NCBI Gene ID: 8569
Uniprot	Q9BUB5
GeneID	8569;
SDS-PAGE MW	51kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HuvEc cells using Mnk1 Polyclonal Antibody



Western Blot analysis of MOUSE-BRAIN cells using Mnk1 Polyclonal Antibody

Note: This product is for in vitro research use only