

Moesin/Ezrin/Radixin Polyclonal Antibody

Catalog No: #41161



Package Size: #41161-1 50ul #41161-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Moesin/Ezrin/Radixin Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	Moesin/Ezrin/Radixin Polyclonal Antibody detects endogenous levels of Moesin/Ezrin/Radixin protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Moesin/Ezrin/Radixin around the non-phosphorylation site of T558.
Target Name	Moesin/Ezrin/Radixin
Other Names	MSN; Moesin; Membrane-organizing extension spike protein; RDX; Radixin; EZR; VIL2; Ezrin; Cytovillin; Villin-2; p81
Accession No.	Swiss-Prot: P26038/P35241/P15311NCBI Gene ID: 4478/5962
Uniprot	P26038
GeneID	4478;
SDS-PAGE MW	67kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of 3T3 cells using Moesin/Ezrin/Radixin Polyclonal Antibody

Note: This product is for in vitro research use only