

## MOR-1 Polyclonal Antibody

Catalog No: #41162

Package Size: #41162-1 50ul #41162-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	MOR-1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	MOR-1 Polyclonal Antibody detects endogenous levels of MOR-1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human MOR-1 around the non-phosphorylation site of S375.
Target Name	MOR-1
Other Names	OPRM1; MOR1; Mu-type opioid receptor; M-OR-1; MOR-1; Mu opiate receptor; Mu opioid receptor; MOP; hMOP
Accession No.	Swiss-Prot: P35372NCBI Gene ID: 4988
Uniprot	P35372
GeneID	4988;
SDS-PAGE MW	60kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/20000.

Not yet tested in other applications.

## Images



Western Blot analysis of 293 cells using MOR-1 Polyclonal Antibody

Note: This product is for in vitro research use only