

MOX-2 Polyclonal Antibody

Catalog No: #41164

Package Size: #41164-1 50ul #41164-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MOX-2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	MOX-2 Polyclonal Antibody detects endogenous levels of MOX-2 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human MOX-2.
Target Name	MOX-2
Other Names	MEOX2; GAX; MOX2; Homeobox protein MOX-2; Growth arrest-specific homeobox; Mesenchyme homeobox 2
Accession No.	Swiss-Prot: P50222NCBI Gene ID: 4223
Uniprot	P50222
GeneID	4223;
SDS-PAGE MW	33kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of COLO205 cells using MOX-2 Polyclonal Antibody



Western Blot analysis of RAT-MUSCLE cells using MOX-2 Polyclonal Antibody

Note: This product is for in vitro research use only