

MRP-L11 Polyclonal Antibody

Catalog No: #41168

Package Size: #41168-1 50ul #41168-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MRP-L11 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF
Species Reactivity	Hu Mk
Specificity	MRP-L11 Polyclonal Antibody detects endogenous levels of MRP-L11 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human MRP-L11.
Target Name	MRP-L11
Other Names	MRPL11; CGI-113; 39S ribosomal protein L11; mitochondrial; L11mt; MRP-L11
Accession No.	Swiss-Prot: Q9Y3B7NCBI Gene ID: 65003
Uniprot	Q9Y3B7
GeneID	65003;
SDS-PAGE MW	21kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

Images



Western Blot analysis of COS7 cells using MRP-L11 Polyclonal Antibody

Note: This product is for in vitro research use only