

MRP-L16 Polyclonal Antibody

Catalog No: #41171

Package Size: #41171-1 50ul #41171-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	MRP-L16 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC
Species Reactivity	Hu
Specificity	MRP-L16 Polyclonal Antibody detects endogenous levels of MRP-L16 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human MRP-L16.
Target Name	MRP-L16
Other Names	MRPL16; PNAS-111; 39S ribosomal protein L16; mitochondrial; L16mt; MRP-L16
Accession No.	Swiss-Prot: Q9NX20NCBI Gene ID: 54948
Uniprot	Q9NX20
GeneID	54948;
SDS-PAGE MW	28kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

Images



Western Blot analysis of 293 cells using MRP-L16 Polyclonal Antibody

Note: This product is for in vitro research use only