

MRP-L24 Polyclonal Antibody

Catalog No: #41174

Package Size: #41174-1 50ul #41174-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MRP-L24 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	MRP-L24 Polyclonal Antibody detects endogenous levels of MRP-L24 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human MRP-L24.
Target Name	MRP-L24
Other Names	MRPL24; 39S ribosomal protein L24; mitochondrial; L24mt; MRP-L24
Accession No.	Swiss-Prot: Q96A35NCBI Gene ID: 79590
Uniprot	Q96A35
GeneID	79590;
SDS-PAGE MW	32kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of HuvEc cells using MRP-L24 Polyclonal Antibody

Note: This product is for in vitro research use only