

MRP-L39 Polyclonal Antibody

Catalog No: #41175

Package Size: #41175-1 50ul #41175-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	MRP-L39 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC
Species Reactivity	Hu
Specificity	MRP-L39 Polyclonal Antibody detects endogenous levels of MRP-L39 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human MRP-L39.
Target Name	MRP-L39
Other Names	MRPL39; C21orf92; MRPL5; RPML5; MSTP003; PRED22; 39S ribosomal protein L39; mitochondrial; L39mt; MRP-L39; 39S ribosomal protein L5, mitochondrial; L5mt; MRP-L5
Accession No.	Swiss-Prot: Q9NYK5NCBI Gene ID: 54148
Uniprot	Q9NYK5
GeneID	54148;
SDS-PAGE MW	39kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

Images



Western Blot analysis of 3T3 cells using MRP-L39 Polyclonal Antibody

Note: This product is for in vitro research use only