

Myf-6 Polyclonal Antibody

Catalog No: #41192

Package Size: #41192-1 50ul #41192-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	Myf-6 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC
Species Reactivity	Hu Ms Rt
Specificity	Myf-6 Polyclonal Antibody detects endogenous levels of Myf-6 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Myf-6.
Target Name	Myf-6
Other Names	MYF6; BHLHC4; MRF4; Myogenic factor 6; Myf-6; Class C basic helix-loop-helix protein 4; bHLHc4; Muscle-specific regulatory factor 4
Accession No.	Swiss-Prot: P23409NCBI Gene ID: 4618
Uniprot	P23409
GeneID	4618;
SDS-PAGE MW	26kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

Images



Western Blot analysis of 293 cells using Myf-6 Polyclonal Antibody

Note: This product is for in vitro research use only