

MYSM1 Polyclonal Antibody

Catalog No: #41197

Package Size: #41197-1 50ul #41197-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MYSM1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu
Specificity	MYSM1 Polyclonal Antibody detects endogenous levels of MYSM1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human MYSM1.
Target Name	MYSM1
Other Names	MYSM1; KIAA1915; Histone H2A deubiquitinase MYSM1; 2A-DUB; Myb-like; SWIRM and MPN domain-containing protein 1
Accession No.	Swiss-Prot: Q5VVJ2NCBI Gene ID: 114803
Uniprot	Q5VVJ2
GeneID	114803;
SDS-PAGE MW	95kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of MCF7 cells using MYSM1 Polyclonal Antibody



Western Blot analysis of MCF7 cells using MYSM1 Polyclonal Antibody

Note: This product is for in vitro research use only