

ND1 Polyclonal Antibody

Catalog No: #41203

Package Size: #41203-1 50ul #41203-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ND1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu
Specificity	ND1 Polyclonal Antibody detects endogenous levels of ND1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human ND1.
Target Name	ND1
Other Names	MT-ND1; MTND1; NADH1; ND1; NADH-ubiquinone oxidoreductase chain 1; NADH dehydrogenase subunit 1
Accession No.	Swiss-Prot: P03886NCBI Gene ID: 4535
Uniprot	P03886
GeneID	4535;
SDS-PAGE MW	36kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of COLO205 cells using ND1 Polyclonal Antibody



Western Blot analysis of L929 cells using ND1 Polyclonal Antibody

Note: This product is for in vitro research use only