

NDUFA4 Polyclonal Antibody

Catalog No: #41205

Package Size: #41205-1 50ul #41205-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NDUFA4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	NDUFA4 Polyclonal Antibody detects endogenous levels of NDUFA4 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human NDUFA4.
Target Name	NDUFA4
Other Names	NDUFA4; NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 4; Complex I-MLRQ; CI-MLRQ; NADH-ubiquinone oxidoreductase MLRQ subunit
Accession No.	Swiss-Prot: O00483NCBI Gene ID: 4697
Uniprot	O00483
GeneID	4697;
SDS-PAGE MW	9kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of RAW264.7 cells using NDUFA4 Polyclonal Antibody

Note: This product is for in vitro research use only