

Neuromedin-S Polyclonal Antibody

Catalog No: #41219



Package Size: #41219-1 50ul #41219-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Neuromedin-S Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC
Species Reactivity	Hu
Specificity	Neuromedin-S Polyclonal Antibody detects endogenous levels of Neuromedin-S protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Neuromedin-S.
Target Name	Neuromedin-S
Other Names	NMS; Neuromedin-S
Accession No.	Swiss-Prot: Q5H8A3NCBI Gene ID: 129521
Uniprot	Q5H8A3
GeneID	129521;
SDS-PAGE MW	24kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.
Immunohistochemistry: 1/100 - 1/300.

Images



Western Blot analysis of MCF7 cells using Neuromedin-S Polyclonal Antibody

Note: This product is for in vitro research use only