

NK-1R Polyclonal Antibody

Catalog No: #41243

Package Size: #41243-1 50ul #41243-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NK-1R Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	NK-1R Polyclonal Antibody detects endogenous levels of NK-1R protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human NK-1R.
Target Name	NK-1R
Other Names	TACR1; NK1R; TAC1R; Substance-P receptor; SPR; NK-1 receptor; NK-1R; Tachykinin receptor 1
Accession No.	Swiss-Prot: P25103NCBI Gene ID: 6869
Uniprot	P25103
GeneID	6869;
SDS-PAGE MW	50kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/5000.

Not yet tested in other applications.

Images



Western Blot analysis of MOUSE-BRAIN COLO cells using NK-1R Polyclonal Antibody

Note: This product is for in vitro research use only