

Nkx-2.6 Polyclonal Antibody

Catalog No: #41244

Package Size: #41244-1 50ul #41244-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nkx-2.6 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	Nkx-2.6 Polyclonal Antibody detects endogenous levels of Nkx-2.6 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Nkx-2.6.
Target Name	Nkx-2.6
Other Names	NKX2-6; NKX2F; Homeobox protein Nkx-2.6; Homeobox protein NK-2 homolog F
Accession No.	Swiss-Prot: A6NCS4NCBI Gene ID: 18092
Uniprot	A6NCS4
GeneID	137814;
SDS-PAGE MW	28kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of COLO205 cells using Nkx-2.6 Polyclonal Antibody

Note: This product is for in vitro research use only