

NM23-H1 Polyclonal Antibody

Catalog No: #41245



Package Size: #41245-1 50ul #41245-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NM23-H1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	NM23-H1 Polyclonal Antibody detects endogenous levels of NM23-H1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human NM23-H1.
Target Name	NM23-H1
Other Names	NME1; NDPKA; NM23; Nucleoside diphosphate kinase A; NDK A; NDP kinase A; Granzyme A-activated DNase; GAAD; Metastasis inhibition factor nm23; Tumor metastatic process-associated protein; nm23-H1
Accession No.	Swiss-Prot: P15531NCBI Gene ID: 4830
Uniprot	P15531
GeneID	4830;
SDS-PAGE MW	23kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HeLa cells using NM23-H1 Polyclonal Antibody



Western Blot analysis of 293T HELA cells using NM23-H1 Polyclonal Antibody

Note: This product is for in vitro research use only