

Nop132 Polyclonal Antibody

Catalog No: #41248

Package Size: #41248-1 50ul #41248-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	Nop132 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC
Species Reactivity	Hu
Specificity	Nop132 Polyclonal Antibody detects endogenous levels of Nop132 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Nop132.
Target Name	Nop132
Other Names	NOL8; C9orf34; NOP132; Nucleolar protein 8; Nucleolar protein Nop132
Accession No.	Swiss-Prot: Q76FK4NCBI Gene ID: 55035
Uniprot	Q76FK4
GeneID	55035;
SDS-PAGE MW	130kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

Images



Western Blot analysis of 3T3 cells using Nop132 Polyclonal Antibody



Western Blot analysis of NIH-3T3 cells using Nop132 Polyclonal Antibody

Note: This product is for in vitro research use only