

Olfactory receptor 8S1 Polyclonal Antibody

Catalog No: #41286

Package Size: #41286-1 50ul #41286-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Olfactory receptor 8S1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	Olfactory receptor 8S1 Polyclonal Antibody detects endogenous levels of Olfactory receptor 8S1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Olfactory receptor 8S1.
Target Name	Olfactory receptor 8S1
Other Names	OR8S1; Olfactory receptor 8S1
Accession No.	Swiss-Prot: Q8NH09NCBI Gene ID: 341568
Uniprot	Q8NH09
GeneID	341568;
SDS-PAGE MW	40kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HT29 cells using Olfactory receptor 8S1 Polyclonal Antibody

Note: This product is for in vitro research use only