

alpha Glucosidase II antibody

Catalog No: #22651



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	alpha Glucosidase II antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC IF
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 47 and 261 of Human GANAB
Target Name	alpha Glucosidase II
Accession No.	Swiss-Prot:Q14697Gene ID:23193
Uniprot	Q14697
GeneID	23193;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 109kd
Western blotting: 1:500-1:3000
Immunohistochemistry: 1:50-1:500
Immunofluorescence: 1:100-1:200


Images



Sample(30 ug whole cell lysate)
A: 293T
B: A431
C: H1299
D: HeLa S3
E: Hep G2
F: MOLT4
G: Raji
7.5% SDS PAGE
Primary antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded SW480 xenograft, using alpha Glucosidase II antibody at 1: 500 dilution.



Confocal immunofluorescence analysis (Olympus FV10i) of methanol-fixed HeLa, using alpha Glucosidase II antibody (Green) at 1: 500 dilution and alpha-tubulin antibody (Red) at 1: 2000.

Note: This product is for in vitro research use only