

PAF65 α Polyclonal Antibody

Catalog No: #41320

Package Size: #41320-1 50ul #41320-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PAF65 α Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms
Specificity	PAF65 α Polyclonal Antibody detects endogenous levels of PAF65 α protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human PAF65 α .
Target Name	PAF65 α
Other Names	TAF6L; PAF65A; TAF6-like RNA polymerase II p300/CBP-associated factor-associated factor 65 kDa subunit 6L; PCAF-associated factor 65-alpha; PAF65-alpha
Accession No.	Swiss-Prot: Q9Y6J9NCBI Gene ID: 10629
Uniprot	Q9Y6J9
GeneID	10629;
SDS-PAGE MW	68kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of 293 cells using PAF65α Polyclonal Antibody



Western Blot analysis of various cells using PAF65α Polyclonal Antibody

Note: This product is for in vitro research use only