

GRSF1 Antibody

Catalog No: #47784

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	GRSF1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	WB, IHC
Species Reactivity	Hu, Ms
Specificity	The antibody detects endogenous levels of total GRSF1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide of human GRSF1
Target Name	GRSF1
Accession No.	Swiss-Prot#:Q12849NCBI Gene ID:2926Gene Accssion:NP_002083
Uniprot	Q12849
GeneID	2926;
Calculated MW	53 kDa
Concentration	0.7
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details


WB dilution:1:500-2000

IHC dilution:1: 30-150


Images



Gel: 8%SDS-PAGE, Lysate: 40 µg, Lane 1-2: 293T and Jurkat cell lysates, Primary antibody:47784(GRSF1 Antibody) at dilution 1/300, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 40 seconds



The image is immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using 47784(GRSF1 Antibody) at dilution 1/25.(Original magnification: 200)



The image is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using 47784(GRSF1 Antibody) at dilution 1/25.(Original magnification: 200)

Background

The protein encoded by this gene is a cellular protein that binds RNAs containing the G-rich element. The protein is localized in the cytoplasm, and has been shown to stimulate translation of viral mRNAs in vitro. Multiple transcript variants encoding different isoforms have been found for this gene.

Note: This product is for in vitro research use only