

NCOA2 (phospho-Ser736) Antibody

Catalog No: #13299

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NCOA2 (phospho-Ser736) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Specificity	Recognizes endogenous levels of NCOA2 (phospho-Ser736) protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human NCOA2.
Target Name	NCOA2
Other Names	BHLHE75; SRC2; TIF2; Nuclear receptor coactivator 2; NCoA-2; Class E basic helix-loop-helix protein 75; bHLHe75; Transcriptional intermediary factor 2; hTIF2
Accession No.	Swiss-Prot#:Q15596NCBI Gene ID:10499
Uniprot	Q15596
GeneID	10499;
Calculated MW	160KD
Concentration	1 mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C


Application Details

Western blotting:1:500 - 1:1000Immunohistochemistry:1:50 - 1:100

Images



Western blot analysis of NCOA2 (phospho-Ser736) expression in HeLa TSA treated whole cell lysates.



Immunohistochemical analysis of NCOA2 (phospho-Ser736) staining in human breast carcinoma formalin fixed paraffin embedded tissue section. The section was then incubated with the antibody at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.

Note: This product is for in vitro research use only