

AKT1/2 (phospho-Thr308/309) Antibody

Catalog No: #13311

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	AKT1/2 (phospho-Thr308/309) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu,Ms,Rt
Specificity	Recognizes endogenous levels of AKT1/2 (phospho-Thr308/309) protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human AKT1/2.
Target Name	AKT1; AKT2
Other Names	PKB; RAC; RAC-alpha serine/threonine-protein kinase; Protein kinase B; PKB; Protein kinase B alpha; PKB alpha; Proto-oncogene c-Akt; RAC-PK-alpha
Accession No.	Swiss-Prot#:P31749; P31751NCBI Gene ID:207; 208
Uniprot	P31749
GeneID	207;
Calculated MW	60KD
Concentration	1 mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting:1:500 - 1:1000

Images



Western blot analysis of AKT1/2 (phospho-Thr308/309) expression in HeLa, PC3 whole cell lysates.

Note: This product is for in vitro research use only