

## GALR2 Antibody

Catalog No: #47876

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	GALR2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu,Ms,Rt
Specificity	Recognizes endogenous levels of GALR2 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human GALR2.
Target Name	GALR2
Other Names	GALNR2; Galanin receptor type 2; GAL2-R; GALR-2
Accession No.	Swiss-Prot#:O43603NCBI Gene ID:8811
Uniprot	O43603
GeneID	8811;
Calculated MW	42KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

## Application Details

WB (1/500 - 1/1000)

## Images



Western blot analysis of GALR2 expression in Hela (A), HepG2 (B), mouse brain (C), mouse kidney (D), rat kidney (E) whole cell lysates

Note: This product is for in vitro research use only