

Ribosomal Protein L10L Polyclonal Antibody

Catalog No: #41399

Package Size: #41399-1 50ul #41399-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ribosomal Protein L10L Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Mk
Specificity	Ribosomal Protein L10L Polyclonal Antibody detects endogenous levels of Ribosomal Protein L10L protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Ribosomal Protein L10L.
Target Name	Ribosomal Protein L10L
Other Names	RPL10L; 60S ribosomal protein L10-like
Accession No.	Swiss-Prot: Q96L21NCBI Gene ID: 140801
Uniprot	Q96L21
GeneID	140801;
SDS-PAGE MW	20kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of HeLa cells using Ribosomal Protein L10L Polyclonal Antibody

Note: This product is for in vitro research use only