

## THNSL2 Antibody

Catalog No: #47957

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	THNSL2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu,Ms,Rt
Specificity	Recognizes endogenous levels of THNSL2 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human THNSL2.
Target Name	THNSL2
Other Names	Threonine synthase-like 2; TSH2; Secreted osteoclastogenic factor of activated T-cells; SOFAT
Accession No.	Swiss-Prot#:Q86YJ6NCBI Gene ID:55258
Uniprot	Q86YJ6
GeneID	55258;
Calculated MW	54KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

## Application Details

WB (1/500 - 1/2000)

## Images



Western blot analysis of THNSL2 expression in HepG2 (A), RAW264.7 (B), H9C2 (C) whole cell lysates

Note: This product is for in vitro research use only